



MEDI CAPTURE

# Medical Imaging Made Easy™



## MediCap™ USB170

The MediCap USB170 is the medical industry's easiest full-featured image recording solution. Your images are saved directly to a USB flash drive or optional hard drive — a much simpler method than trying to capture images with a DVD recorder.



### Medical Image Recorder

The MediCap is even simpler than the pocket camera you use on vacation. It records high quality still images from your endoscope, arthroscope, surgical camera, ultrasound, etc.



### Intuitive Use

It's so easy, most users never look at the manual. Just plug it in and go. You'll be saving images within minutes of opening the box.



### Medical Grade

With its steel case and fluid resistant front panel, the USB170 is built to withstand the rigors of the OR, and it's so compact it will fit on any cart.

# Powerful Options

compatible – versatile – affordable



## USB Hard Drive

Perfect for archiving thousands of images on a single device



## USB Keyboard

Input patient data to organize your image library




## Network Adapter

Turn your USB170 into a mapped drive on your network



## Foot Pedal

Connect your USB170 to your scope's trigger or to a foot pedal



## Also Available...



## MediCap™ USB200 & USB300

Need more than still image recording? Add video recording with the USB200 (standard definition) and the USB300 (high definition).



## MediCap™ USB170 Specs

<b>Device Type</b>	<b>Medical Image Recorder</b> Records still images from standard medical video sources to a USB flash drive or optional USB hard drive
<b>Case</b>	Rugged metal case with high-impact plastic faceplate
<b>LCD Monitor</b>	2 1/2", color (for video confirmation only, not for diagnosis)
<b>Buttons</b>	Sealed membrane, fluid resistant
<b>Video In/Out</b>	PAL or NTSC format using S-video or composite BNC
<b>USB Support</b>	USB2.0, high speed, 1 front connector, 3 rear connectors (all powered)
<b>Networking</b>	Ethernet T10/100 connection via Network Upgrade Kit (optional)
<b>Footswitch/Trigger</b>	Hands-free capture Standard 3.5mm mini jack
<b>Mechanical</b>	Size: 240 x 200 x 63mm / 9.5 x 7.9 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.)
<b>Environmental</b>	Storage: -40° to +85° C Operating: -20° to +40° C
<b>Power</b>	AC 100-240V 50/60Hz 3-prong AC jack
<b>Media</b>	USB flash drives or external USB hard drive
<b>Image Formats</b>	JPEG, TIFF, PNG, DICOM
<b>Image Size</b>	1024 x 768, 800 x 600, 640 x 480 pixels
<b>Certifications</b>	CE 93/42/EEC, UL 60601-1, IEC 60601-1, IEC 60601-1-2, CSA C22.2 No.601.1, E 306417, FCC Part 15, FDA Class 1, RoHS
<b>Kit Contents</b>	MediCap USB170 USB flash drive Video cables Power cable User guide

## MediCapture, Inc.

580 W. Germantown Pike, Ste. 103  
Plymouth Meeting, PA 19462 USA

(888) 922-7887  
(503) 445-6935

www.medicapture.com  
sales@medicapture.com