

MEDI CAPTURE

Medical Imaging Made Easy™



MediCap™ USB200

The MediCap USB200 is the medical industry's easiest full-featured video recording solution. Your videos and still images are recorded directly to a USB flash drive or optional hard drive. The recorded videos are DVD-quality, but without all the hassles of DVD disks.



Medical Video Recorder

Records from your endoscope, arthroscope, surgical camera, ultrasound, etc. Saves to a flash drive or hard drive that can be accessed from any computer on your network.



Intuitive Use

It's so easy, most users never look at the manual. Just plug it in and go. You'll be recording videos within minutes of opening the box.



Medical Grade

With its steel case and fluid resistant front panel, the USB200 is built to withstand the rigors of the OR, and it's so compact it will fit on any cart.

Powerful Options

compatible – versatile – affordable



USB Hard Drive

Perfect for archiving
100+ hours of video
on a single device



USB Keyboard

Input patient data to
organize your photos
and videos





Network Adapter

Turn your USB200
into a mapped drive
on your network



Foot Pedal

Connect your USB200
to your scope's trigger
or to a foot pedal



Also Available...



MediCap™ USB300

Step up to high definition. The USB300 has all the same great features as the USB200, plus HD recording — at a surprisingly affordable price.



MediCap™ USB200 Specs

Device Type	Medical Video Recorder Records video and still images from standard medical video sources to a USB flash drive or optional USB hard drive
Case	Rugged metal case with high-impact plastic faceplate
LCD Monitor	2 1/2", color (for video confirmation only, not for diagnosis)
Speaker	200Hz-15kHz, 1W, (for audio confirmation only, not for diagnosis)
Buttons	Sealed membrane, fluid resistant
Video In/Out	PAL or NTSC format using S-video or composite BNC
Audio In/Out	Standard 3.5mm stereo phone jacks
USB Support	USB2.0, high speed, 1 front connector, 3 rear connectors (all powered)
Networking	Ethernet T10/100 connection via Network Upgrade Kit (optional)
Footswitch/Trigger	Hands-free capture Standard 3.5mm mini jack
Mechanical	Size: 240 x 200 x 63mm / 9.5 x 7.9 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.)
Environmental	Storage: -40° to +85° C Operating: -20° to +40° C
Power	AC 100-240V 50/60Hz 3-prong AC jack
Media	USB flash drives or external USB hard drive
Image Formats	JPEG, TIFF, PNG, DICOM
Image Size	1024 x 768, 800 x 600, 640 x 480 pixels
Video Format	MPEG2 (DVD-quality)
Audio Format	MPEG1 L2, 2 channels
Certifications	CE 93/42/EEC, UL 60601-1, IEC 60601-1, IEC 60601-1-2, CSA C22.2 No.601.1, E 306417, FCC Part 15, FDA Class 1, RoHS
Kit Contents	MediCap USB200 USB flash drive Video cables Power cable User guide

MediCapture, Inc.

580 W. Germantown Pike, Ste. 103
Plymouth Meeting, PA 19462 USA

(888) 922-7887
(503) 445-6935

www.medicapture.com
sales@medicapture.com